

Circular 006/2025

20th February 2025

Dear Resident

MUIZENBERG - TELECARE Medical Emergency Procedure

Telecare is our medical response service provider and in order to ensure that we are well prepared in the event of a medical emergency, I thought a reminder of our emergency protocols would be in order.

TELECARE UNIT

ALL homes are equipped with a Telecare unit which is linked to a 24/7 call centre. This is the device to use in any medical emergency not the internal phone. You will have been issued with a panic button and strongly urge you to wear this wrist band or necklace on your person at all times.

In a medical emergency you should firmly <u>press the red button on your Telecare unit or</u> <u>the red button on your Telecare panic button</u> and the operator should answer your call within 60 seconds. They will establish your emergency. As they already have your medical history, they will call the paramedics and ambulance, if required and brief them en-route.


Once the Telecare call centre has established that there is a medical emergency, they will also call the on-duty Duty Manager to go to your unit and assist and wait with you for the ambulance.

Do not use your landline whilst waiting for a response from the Telecare Centre. Only press the red button **ONCE** and do not press any other button on your Telecare unit, <u>as it cancels</u> the emergency call.

Please note that the health staff and private carers on site (other than the Sisters or RN's on duty) may not be equipped to address ALL emergencies and this may preclude them from performing certain medical procedures. They are trained in basic life support which may include first aid and CPR. They have been given the mandate to intervene or commence such treatment until the appropriate emergency service provider arrives and resumes care.

While our staff check the functioning of Telecare units/internal phone connections routinely, it's wise to test it regularly yourself. We encourage you to test the system from time-to-time by pushing the red button. This is not "crying wolf", it is an important safety check on the equipment.

Simply tell the Telecare operator this when they respond - "Thank you, I am just testing my Telecare response/unit". Report any failed test to reception for immediate attention. I have attached a Telecare "how to" diagram. Please request a copy from reception.



We would like to request that you wear your blue buddy band, especially if you are going out of the village to the shops. The buddy band has a unique number which identifies that you are a resident at Evergreen Muizenberg, it has our contact details as well as any medical conditions and lists your next-of-kin.

Do not hesitate to call Telecare because you think that the emergency is not severe enough. No harm is done if your condition turns out to be a false alarm.

There is more risk in not calling them and losing critical time. Please feel free to call me if you need any additional information or wish to discuss any aspect of the emergency protocol.

Yours sincerely

Christine Dempers Village Manager – Muizenberg

^{\\}amdeccpt\evergreen\Muizenberg\1-MEMO & NEWSLETTERS (Communications)\2025\20250127 - Circular 006 - Telecare.Emergency Response Update.docx